

Basic Steam Sterilizers and Specialized Stainless Steel Products

Main Parts Utilized

Serve Health Forever



Various Specialized Stainless Steel Products



CSSD Aired Shelving

Vented open-wire design	1200x450x180mm(LxWxH)
Minimum Condensate Formation	1200x600x180mm(LxWxH)
Extended Shelf Life	900x450x180mm(LxWxH)
High Tec Design	900x600x180mm(LxWxH)



Horizontal round pressure steam sterilizer

Sterilizing Chamber Size	Φ 600x1000mm
Physical Size	1360x760x1820mm(LxWxH)
Working Pressure	0.145MPa~0.165MPa
Working Temperature	126 °C~129 °C
Heat Equilibrium	<±1°C
Power Source	12KW/AC380V.50Hz
Steam Source	0.3-0.8MPa
Net Weight: Steam Heating	335Kg
Electric Heating Type	400Kg
Electric/Steam Heating	410Kg

Leadership Trough Technology



Horizontal rectangle pressure steam sterilizer

Chamber Size Single Door	1420x760x1065mm(LxWxH)
Chamber Size Double Door	1520x760x1065mm(LxWxH)
Physical Size Single Door	1900x1200x2000mm(LxWxH)
Physical Size Double Door	2100x1200x2000mm(LxWxH)
Working Pressure	0.145MPa
Working Temperature	126 °C
Heat Equilibrium	<±1°C
Steam Source	0.3-0.8MPa
Net Weight: Single Door	1600Kg
Net Weight: Double Door	1860Kg

